

GUS Histochemical Staining (w/ Acetone fixation)

2025-06-17 Woo-Taek Jeon

Reagent / condition	Preparation or use
50 mL, 3 mM ferri/ferro	46.5 mL base solution + 1.5 mL ferri (3 mM) + 1.5 mL ferro (3 mM) + 0.5 mL X-Gluc (1 mM)
50 mL, 10 mM ferri/ferro	39.5 mL base solution + 5 mL ferri (10 mM) + 5 mL ferro (10 mM) + 0.5 mL X-Gluc (1 mM)

1. Fix tissue in ice-cold 90% acetone for 2 h.
2. Alternatively, infiltrate under vacuum at room temperature for 5 min, release the vacuum and incubate for 20-30 min at room temperature, then return samples to ice.
3. Infiltrate samples in staining solution without X-Gluc for 5 min under vacuum on ice.
4. Replace with X-Gluc staining solution and infiltrate for 5 min under vacuum on ice in the dark.
5. Incubate at 37 °C in the dark until the desired signal develops.

Source: [GUS reference](#)